

AB Sample Kit

Urine Sample Splitting Kit

QUALITY
ISO 9001

ENVIRONMENT
ISO 14001

HEALTH & SAFETY
ISO 45001

A-B Sample Kit for Splitting Urine Samples

The Australian Standard for urine testing (AS4308) dictates that you must split the sample in the event you get an unconfirmed positive at screening.

We have created a simple collection kit making the splitting process easy and clean.

To split the sample:

- Simply insert transfer tube into sample,
- Insert vacuette sample tube into transfer tube.
- Sign, seal and send.

A-B Split Sample Kit Features:

- Vacuette tubes to automatically suck up sample into collection tube
- 3 Sets of serial numbered security seals
- Biohazard transport bag
- Absorbant Strip



KEY DETAILS

- Meets AS4308 requirements
- No pouring urine samples
- Quality Guaranteed
- Simple process
- Economical

Call Today!
1300 79 70 30