

# SureStep 7 + 6

Multi-Drug Test panel with PINACA Test (K3 Syncann)

**QUALITY**  
ISO 9001

**ENVIRONMENT**  
ISO 14001

**HEALTH & SAFETY**  
ISO 45001



## Detect \*PINACA\* with Ease

The SureStep urine drug screening device provides a quick and simple mechanism to test for multiple drugs of abuse with the added benefit of a built in panel for checking sample validity.

The device is the latest in the Surestep device stable and detects the latest synthetic cannabis (Syncann) group (PINACA or K3).

This group of syncanns are very potent and dangerous in the workplace.

Detection with this device is easy and fast.

With negative results indicated in a matter of minutes and detected results at 5 minutes this is your industry-best onsite test device for minimal downtime and effective detection rates.

## KEY DETAILS

- Onsite Screening Device
- Designed for Workplace
- Tests for PINACA
- Accurate & Reliable
- Proven 3 Step Process
- Tests for 7 Drugs of Abuse & 6 Adulterants

**Call Today!**  
**1300 79 70 30**


SureStep 7 + 6 Drug Test Cup – Urine instant drug test with built in adulteration test panel provides a rigorous and affordable workplace drug test solution that is aligned to Australian Standards.

The SureStep 7 + 6 Drug Test Cup offers a unique advanced split key design that allows the test operator to activate the test and split 3.5mL of urine from the main volume of sample and into a test chamber – all at the push of a button.

No handling of urine required in order to complete the test and the integrity of the sample collected is not compromised and is available for further analysis at a laboratory if required.

An inbuilt thermometer to check and verify the temperature of the sample combined with an adulteration test panel to check the integrity are key features of this smart advanced technology.

CMM Technology™ offer free online training: we have developed an online product certification program for our clients to use with this product so that you can rest assured your staff are fully conversant in the use of this device.



## Key Product Information:

- One-step – no urine handling
- Tests for the AS4308 drugs of abuse & PINACA
  - **Amphetamines**
  - **Benzodiazepines**
  - **Cocaine**
  - **Methamphetamine**
  - **Marijuana (THC)**
  - **PINACA (K3)**
  - **Opiates**
- Adulteration Test Panel
  - **Creatinine**
  - **Nitrite**
  - **Gluteraldehyde**
  - **pH**
  - **Specific gravity**
  - **Oxidants/PCE**
- Test starts when the collector is ready (not when the subject gives the specimen)
- Built-in thermometer
- Easy open/close no-leak container
- Only 25ml of urine is needed (NB: for sample confirmation, only 3.5ml of urine is required to actually run the test)
  - Splits the instant test sample from the main laboratory sample
- Read results within 2-5 minutes
- Kits can be tailored to suit your specific requirements