

## LATEST TECHNOLOGY

ToxWipe 7 Oral Fluid Drug Screen Device is developed to detect AS4760:2019 listed active drugs-of abuse present in saliva.

The simple collection procedure makes oral fluid testing a preferred method for many.

ToxWipe 7 Features:

- Performance independently verified by an AS4760:2019 accredited laboratory.
- No handling or collection of saliva.
- Utilizes a very small volume to activate and test.
- Integrated buffer.
- Patented collection pad.
- Latest generation testing antibodies.

Tests for the following AS4760:2019 Drug Classes:

DRUG CLASS	CUT-OFF
Amphetamine	50
Benzodiazepines	50
Cocaine	50
Methamphetamine	50
Opiates	50
Oxycodone	40
THC (Cannabis)	15



## KEY DETAILS

- Results in minutes
- Tests for 7 drugs of abuse
- Quality Guaranteed
- AS4760:2019 Drugs
- Accurate
- Reliable
- Proven
- No oral fluid handling

**Call Today!**  
**1300 79 70 30**